



Ref.No.1-23/Pur/Web/ICMR-RMRCSP/ 386

Date: 06-February-2026

खरीद सूचना/PROCUREMENT NOTICE

ICMR-रीजनल मेडिकल रिसर्च सेंटर, श्री विजया पुरम, प्रतिष्ठित, योग्य और अधिकृत सप्लायर/विक्रेताओं से उन वस्तुओं की सप्लाई और डिलीवरी के लिए कोटेशन आमंत्रित करता है, जो GeM प्लेटफॉर्म पर उपलब्ध नहीं हैं, सेंटर के कॉम्प्लेक्स, डॉलीगंज, श्री विजया पुरम में, जैसा कि नीचे बताया गया है:

ICMR-Regional Medical Research Centre, Sri Vijaya Puram invites quotations from reputed, eligible, and authorized suppliers/vendors for supply and delivery of the items, which are not available in GeM platform, at the centre's complex at Dollygunj, Sri Vijaya Puram as detailed below:

I. Reference No. - 264191

Serial No.	Name of the items	Quantity (Units)
1.	TAE buffer (50x) 1000 ml	1
2.	Formaldehyde 500 ml	1
3.	MgCl2 (50x) 1ml	5

II. Reference No. - 260509

Sl. No.	Item	Quantity
1	BLUE PEN	1 No.s
2	BLACK PEN	1 No.s
3	RED PEN	1 No.s
4	BLUE MARKER	1 No.s
5	BLACK MARKER	1 No.s
6	CELLOTAPE	5 No.s
7	SCALE	1 No.s
8	CARDBOARD	5 No.s
9	TISSUE	10 No.s
10	BASKET WITH HANDLES	2 No.s

III. Reference No. - 264182

Sl. No.	Item	Quantity	
1	Plastic Container	30 No.s	Plastic Containers with Wide Mouth, Chemical Reagent Storage, transparent, 300-ML
2	Silica gel (500 gm)	1 No.s	white silica gel in crystalline form
3	Storage Basket (Field Containers)	2 No.s	Storage Basket with Handles, Solid Plastic Big Bin, Handles Plus Lid, Size 30 x 41 x 28, Rectangular.

IV. Reference No. – 259850

Serial No.	Name of the items	Quantity
1.	Celling fan	1 Unit
2.	Choke	3 Units
3.	LED slim surface panel light (18 watt make Crompton)	5 Units
4.	LED bulb (18 watt)	5 Units
5.	LED slim surface panel light (8 watt)	5 Units

V. Reference No. - 262317

Serial No.	Name of the items	Quantity (Units)
1.	Panel light, 2*2 feet made by Crompton/Panasonic	15
2.	Panel light, Round,20 watt, made by Crompton/Panasonic	20
3.	Fan, ceiling (bajaj) white	20
4.	Mirror light ,2 feet made by Crompton/Panasonic	10
5.	Down light, square,8W, made by Crompton/Panasonic	30
6.	Down light, Round,18W, made by Crompton/Panasonic	10
7.	Down light, sound, 2*2 feet, made by Crompton/Panasonic	05
8.	Down light, Round,8W, made by Crompton/Panasonic	10
9.	LED light,15W, made by Panasonic/Bajaj	15
10.	Streetlight, LED,35W, made by Crompton/Panasonic	15
11.	Gate light, Big size, made by Crompton/Bajaj	04
12.	Wire 1.5mm, made by Anchor/Havells	01 coil
13.	Insulation tape	10 Rolls

VI. Reference No. - 258034

Serial No.	Name of the items	Quantity
1	Autoclave bag Medium (Transparent)	100 units
2	Autoclave bag Small (Transparent)	100 units
3	Cryolabels (Small)	1000 units

✓ 02/07/2026

4.	Cryolabels (Medium)	1000 units
5.	Glass trough	1 Unit
6.	Plastic basket	3 Units

VII. Reference No. – 264931

Sl. No.	Item	Quantity
1	CFL Lamp ,12 W	8 No.s
2	LED, Panel Light Round 18 W	8 No.s
3	LED, Surface Light Round 24 W	8No.s
4	LED Bulb 18 W	5No.s
5	Celling Fan, White	5No.s
6	Mirror Light, LED 2ft	9No.s
7	LED, Surface light 1x1 ft 36 W	4No.s
8	Bed Lamp	6No.s
9	Tube light Led 36/40W	3No.s
10	Wall Decorative Light	4No.s
11	1.5 Sq.mm copper wire	1 coil
12	Insulation Tape	10No.s

VIII. Reference No- 255309

Sl. No	Name of the Item	Quantity
01	RSV IgG Enzyme Immunoassay Kit	4
02	RSV IgM Enzyme Immunoassay Kit	2
03	RSV IgA Enzyme Immunoassay Kit	2

IX - Reference No. 264190

Sl. No.	Name of the Item	Quantity
1	Micro Pestle (Polypropylene, 82 mm length, tissue homogenizer with beveled tip; suitable for DNA/RNA work	50 os

X. Reference No. – 264183

Serial No.	Name of the items	Quantity
1.	Stickers or Labels	20 sheet

✓ 26/10/2016

	(Size 25 x 12 mm)	
2.	Cello brown tape 3inch	5 Units
3.	Cello white tape 1 inch	45 Units

XI. Reference No. – 264183

Serial No.	Name of the items	Quantity
1.	Dipper or Laddle	4 units
2.	Aspirator (sucking tube)	6 Units
3.	Mosquito Rearing Cages	5 Units
4	Culture tray	10 Units
5	Entomological Forceps	04 units
6	Entomological dissection box	1 Box
7	Water Filter/Strainer	04 units
8	Insect Dissection Needle	2 units
9	Enamel Bowel	04 units

Specifications for items in Reference No. – 264183

Dipper or Laddle:

Fixed-length handle: 60–90 cm (default: 60 cm), Three-segment telescopic aluminium handle, Length option 1: Maximum of 200 cm., Collection bowel should be Enamel material 350 ml, Three-segment telescopic aluminium handle

Aspirator (sucking tube)

L30 x D12mm Acrylic Straight for insectary and field collection of mosquitoes

Mosquito Rearing Cages

W24.5 x D24.5 x H24.5 cm, Netting Nylon Type, Dimensions: 32.5 x 32.5 x 32.5 cm, Lightweight aluminium frame, Fine nylon mesh for ventilation, Easy assembly with plastic connectors, Front sleeve opening for easy access, Solid floor for added stability

Culture tray

Dimensions: L37.5 x W30 x H7.5 cm, Net Weight: 460 grams, Main Material: Plastic- Polypropylene, Mesh Size: 26 x 24 | 1, 220-micron aperture

Entomological Forceps

Feather forceps, L12 cm, 3 grams 9.

W. H. M. 10/12/2022

Insect Dissection Needle - Insect Dissection Needle LI-IR-55 for Entomological work, The sharp tips with varying lengths (5cm or 2.5cm)

Enamel Bowel

16cm (6 inches) in diameter, ideal for medical, laboratory, and clinical use

Interested firms/vendors who can supply the above-mentioned item(s) are requested to submit their quotation (including GST, delivery charges, and warranty details, if applicable) through email.

XII. Reference No- 255312

Sl No	Name of the Item	Quantity
01	Unmodified Oligos of various lengths, (Primers for PCR)	1
02	Dual-labelled Probes of various lengths, TaqMan probes for real-time PCR	1

Specification: Reference No- 255312

S. No.	Oligo name	Sequence (5' to 3')	No. of bases	Synthesis Scale (μmole)	Purification
1	RSV A FP	AGATCAACTTCTGTCATCCAGCAA	24	0.2	DST
2	RSV A RP	TTCTGCACATCATATAATTAGGAG	22	0.2	DST
3	RSV B FP	AAGATGCAAATCATAAATTACAGGA	26	0.2	DST
4	RSV B RP	TGATATCCAGCATCTTAAGTA	22	0.2	DST

Table 1: Detailed Specification for Oligos

Table 2: Detailed Specification for Probes

S. No.	Oligo name	5' Mod	Sequence (5' to 3')	3' Mod	No. of bases	Synthesis Scale (nmole)
1	RSV A P	FAM	CACCATCCAACGGAGCACAGGAGAT	MGBNFQ	25	50
2	RSV B P	VIC	TTTCCCTTCCTAACCTGGACATA	MGBNFQ	23	50

✓ 06/01/2026

XIII. Reference No- 255312

रम सां/ S.No.	कतया/ सतमतन कत पूरा ब्योरत Name of item(s)/work with full specification	सां/Qty
1	Dengue Serotyping-Chikungunya-Zika Multiplex Probe PCR Kit (100R)	1 No
2	cDNA Synthesis Kit (100R)	1 No
3	Viral RNA Kit for RNA Extraction (250R)	1 No

Specification: for Reference No- 255312

Sl. No.	Name of Item	Specification
1	Dengue Serotyping-Chikungunya-Zika Multiplex Probe PCR Kit (100R)	<ul style="list-style-type: none"> Detection: Simultaneous serotyping of Dengue viruses and detection of Chikungunya and Zika in a single assay. Time consumption: Results within 90 minutes. Assay: One-step assay i.e. reverse transcription and amplification are performed in same tube. Reaction: Two Tube reaction to detect DENV1, DENV2, DENV3, DENV4, CHIKV and ZIKA. Product contents: Should include all reagents & controls for validity of the test. Limit of Detection: 2.87 copies/µL for DENV-1, < 1 copy/µL for DENV-2, 1.87 copies/µL for DENV-3, 1.41 copies/µL for DENV-4, < 1 copy/µL for CHIKV and < 1 copy/µL for ZIKV. <p>Compatibility: With 4-channel and 5-channel qPCR cyclers.</p>
2	cDNA Synthesis Kit (100R)	<ul style="list-style-type: none"> Enzyme System: M-MuLV RT: Efficient RNA-dependent DNA polymerase for cDNA synthesis. RNase Inhibitor: Integrated protection against RNA degradation during the reaction. <p>PRIMING OPTIONS</p> <ul style="list-style-type: none"> Random Hexamers: For non-specific priming of all RNA species. Oligo (dT): For specific targeting of mRNA poly(A) tails. Primer Mix: Pre-mixed hexamers and oligo (dT) for balanced library prep. <p>OPERATING CONDITIONS</p> <ul style="list-style-type: none"> Temperature: Optimal synthesis at 42°C; avoid exceeding 50°C and Storage: Stable for 12–15 months when stored at -20°C. <p>PERFORMANCE STATISTICS</p> <ul style="list-style-type: none"> Conversion Rate: High-yield transformation of mRNA into single-stranded cDNA. Buffer System: Optimized for consistent results across varying

26/01/2016

		<p>template concentrations.</p> <p>SPECIFICITY</p> <ul style="list-style-type: none"> • Priming Accuracy: Optimized buffer-enzyme mix reduces non-specific binding. • Target Sensitivity: Enhanced detection when paired with gene-specific primers. <p>QUALITY CONTROL</p> <ul style="list-style-type: none"> • Nuclease Free: Certified free of endonucleases, exonucleases, and RNases. <p>Reproducibility: Batch-tested to ensure single, bright bands in downstream PCR and should be 100%.</p>
3	Viral RNA Kits for RNA Extraction (250R)	<ul style="list-style-type: none"> • Rapid isolation of high-quality, ready-to-use RNA. <p>No organic extraction or alcohol precipitation should be used.</p> <ul style="list-style-type: none"> • Recovery: > 90% (Highly efficient) of Viral RNA at any titre. • Sample type: Cell free supernatant. • Processing: Manual. • Elution volume: 80-100 μl. • Format: Spin columns. • Sample amount: 140 μl. • Time per run or prep: 20-40 minutes. • Removal of complete contaminants and inhibitors. <p>Carrier RNA based extraction for enhanced binding of viral nucleic acids.</p>

Quotation should be sent to : **procurement1rmrcsvp@gmail.com**

Due date and time for submission of quotation: **13 February 2026 by 3.30 pm**

नियम और शर्तें:

1. कोटेशन फर्म के लेटरहेड पर डायरेक्टर ICMR-RMRC श्री विजया पुरम डॉलीगंज अंडमान और निकोबार द्वीप समूह भारत पिनकोड 744103 के नाम से जमा किया जाना चाहिए।
2. कोटेशन ऊपर बताए गए रेफरेंस और सीरियल डिटेल के अनुसार तैयार किया जाना चाहिए।
3. GST नंबर और कोटेशन की वैलिडिटी साफ तौर पर बताई जानी चाहिए।
4. डिलीवरी की अवधि कोटेशन में बताई जानी चाहिए।
5. संस्थान बिना कोई कारण बताए किसी भी या सभी कोटेशन को स्वीकार या अस्वीकार करने का अधिकार सुरक्षित रखता है।
6. तय तारीख और समय के बाद मिले कोटेशन पर विचार नहीं किया जाएगा।
7. कोटेशन और सभी सपोर्टिंग डॉक्यूमेंट/सर्टिफिकेट सिफ़र पासवर्ड प्रोटेक्टेड PDF फ़ॉर्मेट में ही सबमिट किए जाएंगे, दूसरे फ़ाइल फ़ॉर्मेट पर विचार नहीं किया जाएगा।
8. PDF फ़ाइल खोलने का पासवर्ड कृपया 13 फ़रवरी 2026, शाम 3.30 बजे के बाद अलग से ईमेल द्वारा भेजें।
7. दरें ICMR-RMRC CSVP पर सभी टैक्स और डिलीवरी चार्ज सहित होनी चाहिए।

Wifn 09/01/2026

Terms & Conditions:

1. Quotation should be submitted on the firm's letterhead in the name of Director ICMR-RMRC Sri Vijaya Puram, Dollygunj, Andaman & Nicobar Islands India, Pincode:744103
2. Quotation must be prepared as per the Reference and serial detailed above.
3. GST number and validity of quotation must be clearly mentioned.
4. Delivery period should be mentioned in the quotation
5. The institute reserves the right to accept or reject any or all quotations without assigning any reason.
6. Quotations received after the Due date and time will not be considered.
7. The quotation along with all supporting documents/certificates shall only be submitted in **Password protected PDF format**, other file formats will not be considered.
8. The password for opening the PDF file(s) may kindly be sent by email separately after **13 February 2026, 3.30 pm**
9. Rates must be inclusive of all taxes and delivery charges at ICMR-RMRC SVP.

२५/०२/२०२६
प्रशासनिक अधिकारी
Administrative Officer
कृते निदेशक/for Director